

For Research Use Only

GATM Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68513-6-PBS



Basic Information

Catalog Number:

68513-6-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG3754

GenBank Accession Number:

BC004141

GeneID (NCBI):

2628

UNIPROT ID:

P50440

Full Name:

glycine amidinotransferase (L-arginine:glycine amidinotransferase)

Calculated MW:

423 aa, 48 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3E8G1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50565-3, GATM Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68513-4-PBS. Detection antibody: 68513-6-PBS. Standard: Ag3754. Range: 0.781-12.5 ng/mL.