

For Research Use Only

DSG2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68515-2-PBS



Basic Information

Catalog Number: 68515-2-PBS	GenBank Accession Number: BC099655	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 1829	CloneNo.: 1C7E6
Source: Mouse	UNIPROT ID: Q14126	
Isotype: IgG1	Full Name: desmoglein 2	
Immunogen Catalog Number: AG20633	Calculated MW: 1118 aa, 122 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51241-1, DSG2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68515-1-PBS. Detection antibody: 68515-2-PBS. Standard:Ag20633. Range: 0.391-100 ng/mL.