

For Research Use Only

ATL3 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68516-5-PBS



Basic Information

Catalog Number: 68516-5-PBS	GenBank Accession Number: BC077727	Purification Method: Protein G Magarose purification
Concentration: 1 mg/ml	GeneID (NCBI): 25923	CloneNo.: 3D4E10
Source: Mouse	UNIPROT ID: Q6DD88	
Isotype: IgG1	Full Name: atlastin GTPase 3	
Immunogen Catalog Number: AG10248	Calculated MW: 541 aa, 61 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50703-3, ATL3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68516-5-PBS. Detection antibody: 68516-6-PBS. Standard: Ag10248. Range: 0.781-100 ng/mL.