

For Research Use Only

MTHFD2 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68524-2-PBS

Basic Information

Catalog Number:

68524-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG33605

GenBank Accession Number:

BC017054

GeneID (NCBI):

10797

UNIPROT ID:

P13995

Full Name:

methylenetetrahydrofolate
dehydrogenase (NADP+ dependent) 2,
methenyltetrahydrofolate
cyclohydrolase

Calculated MW:

350 aa, 38 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1B7H2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50286-1, MTHFD2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68524-2-PBS. Detection antibody: 68524-3-PBS. Standard: Ag33605. Range: 0.391-100 ng/mL.