

For Research Use Only

UBA5 Monoclonal antibody, PBS Only



Catalog Number: 68537-1-PBS

Basic Information

Catalog Number:

68537-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG33685

GenBank Accession Number:

BC009737

GeneID (NCBI):

79876

UNIPROT ID:

Q9GZZ9

Full Name:

ubiquitin-like modifier activating enzyme 5

Calculated MW:

404 aa, 45 kDa

Observed MW:

45 kDa

Purification Method:

Protein A purification

CloneNo.:

3E5C8

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human, mouse, rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68537-1-Ig (UBA5 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68537-1-PBS in a different storage buffer formulation.