

For Research Use Only

HLA-DRA Monoclonal antibody, PBS Only



Catalog Number: 68543-1-PBS

Basic Information

Catalog Number: 68543-1-PBS	GenBank Accession Number: BC032350	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 3122	CloneNo.: 3D6E8
Source: Mouse	ENSEMBL Gene ID: ENSG00000204287	
Isotype: IgG1	UNIPROT ID: P01903	
Immunogen Catalog Number: AG11224	Full Name: major histocompatibility complex, class II, DR alpha	
	Calculated MW: 254 aa, 29 kDa	
	Observed MW: 32-34 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human, Pig

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68543-1-Ig (HLA-DRA antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68543-1-PBS in a different storage buffer formulation.