

For Research Use Only

KRT12 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68545-3-PBS



Basic Information

Catalog Number: 68545-3-PBS	GenBank Accession Number: BC156641	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 3859	CloneNo.: 2E9A8
Source: Mouse	UNIPROT ID: Q99456	
Isotype: IgG1	Full Name: keratin 12	
Immunogen Catalog Number: AG20332	Calculated MW: 494 aa, 54 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51143-1, KRT12 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68545-3-PBS. Detection antibody: 68545-4-PBS. Standard:Ag20332. Range: 0.098-12.5 ng/mL



Cytometric bead array standard curve of MP51143-3, KRT12 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68545-3-PBS. Detection antibody: 68545-7-PBS. Standard:Ag20332. Range: 0.098-1.563 ng/mL