

For Research Use Only

# SUCLA2 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68549-2-PBS

## Basic Information

<b>Catalog Number:</b> 68549-2-PBS	<b>GenBank Accession Number:</b> BC027587	<b>Purification Method:</b> Protein G purification
<b>Size:</b> 1mg/ml	<b>GeneID (NCBI):</b> 8803	<b>CloneNo.:</b> 2C11G2
<b>Source:</b> Mouse	<b>UNIPROT ID:</b> Q9P2R7	
<b>Isotype:</b> IgG1	<b>Full Name:</b> succinate-CoA ligase, ADP-forming, beta subunit	
<b>Immunogen Catalog Number:</b> AG17967	<b>Calculated MW:</b> 463 aa, 50 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50128-1, SUCLA2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68549-2-PBS. Detection antibody: 68549-3-PBS. Standard: Ag17967. Range: 3.125-100 ng/mL.