

For Research Use Only

RAB35 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68552-4-PBS



Basic Information

Catalog Number: 68552-4-PBS	GenBank Accession Number: BC015931	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 11021	CloneNo.: 2G8E1
Source: Mouse	UNIPROT ID: Q15286	
Isotype: IgG2a	Full Name: RAB35, member RAS oncogene family	
Immunogen Catalog Number: AG21424	Calculated MW: 201 aa, 23 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50998-2, RAB35 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68552-4-PBS. Detection antibody: 68552-3-PBS. Standard:Ag21424. Range: 0.391-100 ng/mL.