

For Research Use Only

# GSPT1 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68553-2-PBS

## Basic Information

Catalog Number: 68553-2-PBS	GenBank Accession Number: BC009503	Purification Method: Protein G purification
Size: 1mg/ml	GeneID (NCBI): 2935	CloneNo.: 1H5E8
Source: Mouse	UNIPROT ID: P15170	
Isotype: IgG1	Full Name: G1 to S phase transition 1	
Immunogen Catalog Number: AG27986	Calculated MW: 68 aa, 4 kDa	

## Applications

Tested Applications:  
Sandwich ELISA, Indirect ELISA

Species Specificity:  
Human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sandwich ELISA standard curve of MP50094-1, GSPT1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68553-2-PBS. Detection antibody: HRP-conjugated 68553-3-PBS. Standard: Ag27986. Range: 0.078-10 ng/mL.