

For Research Use Only

ALAS1 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68567-3-PBS

Basic Information

Catalog Number:

68567-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG33580

GenBank Accession Number:

BC011798

GeneID (NCBI):

211

UNIPROT ID:

P13196

Full Name:

aminolevulinate, delta-, synthase 1

Calculated MW:

640 aa, 71 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2E1E3

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50232-1, ALAS1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68567-2-PBS. Detection antibody: 68567-3-PBS. Standard: Ag33580. Range: 0.195-25 ng/mL