

For Research Use Only

# PIP5K3 Monoclonal antibody, PBS Only



Catalog Number: 68571-1-PBS

## Basic Information

Catalog Number:

68571-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG4366

GenBank Accession Number:

BC032389

GeneID (NCBI):

200576

UNIPROT ID:

Q9Y2I7

Full Name:

phosphatidylinositol-3-phosphate/phosphatidylinositol 5-kinase, type III

Calculated MW:

451 aa, 50 kDa, 237 kDa

Observed MW:

240-250 kDa

Purification Method:

Protein G purification

CloneNo.:

2F5G1

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human, mouse, rat, rabbit, pig

## Background Information

## Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68571-1-Ig (PIP5K3 antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68571-1-PBS in a different storage buffer formulation.



Various lysates were subjected to SDS PAGE followed by western blot with 68571-1-Ig (PIP5K3 antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68571-1-PBS in a different storage buffer formulation.