

For Research Use Only

CYP46A1 Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 68576-4-PBS



Basic Information

Catalog Number:

68576-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG16786

GenBank Accession Number:

BC022539

GeneID (NCBI):

10858

UNIPROT ID:

Q9Y6A2

Full Name:

cytochrome P450, family 46, subfamily A, polypeptide 1

Calculated MW:

500 aa, 57 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1G10C8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50891-2, CYP46A1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68576-4-PBS. Detection antibody: 68576-5-PBS. Standard:Ag16786. Range: 0.391-100 ng/mL



Cytometric bead array standard curve of MP50891-3, CYP46A1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68576-6-PBS. Detection antibody: 68576-4-PBS. Standard:Ag16786. Range: 0.391-100 ng/mL