

For Research Use Only

PTPIP5 1 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68582-3-PBS

Basic Information

Catalog Number: 68582-3-PBS	GenBank Accession Number: BC008970	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 55177	CloneNo.: 1F3H2
Source: Mouse	UNIPROT ID: Q96TC7	
Isotype: IgG1	Full Name: family with sequence similarity 82, member A2	
Immunogen Catalog Number: AG15135	Calculated MW: 470 aa, 52 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50287-1, PTPIP51 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68582-2-PBS. Detection antibody: 68582-3-PBS. Standard: Ag15135. Range: 0.781-50 ng/mL