

For Research Use Only

SSH3 Monoclonal antibody

Catalog Number: 68583-1-Ig



Basic Information

Catalog Number:

68583-1-Ig

Size:

1000 µg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG13139

GenBank Accession Number:

BC004210

GeneID (NCBI):

54961

UNIPROT ID:

Q8TE77

Full Name:

slingshot homolog 3 (Drosophila)

Calculated MW:

659 aa, 73 kDa

Observed MW:

90-95 kDa

Purification Method:

Protein A purification

CloneNo.:

1A9C9

Recommended Dilutions:

WB 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Positive Controls:

WB : A431 cells, A549 cells, T-47D cells, MCF-7 cells, HeLa cells, HT-29 cells, HepG2 cells, HuH-7 cells, THP-1 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68583-1-Ig (SSH3 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.