

For Research Use Only

BCAT2 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68586-4-PBS

Basic Information

Catalog Number:

68586-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG9595

GenBank Accession Number:

BC001900

GeneID (NCBI):

587

UNIPROT ID:

O15382

Full Name:

branched chain aminotransferase 2, mitochondrial

Calculated MW:

44 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3H1B3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50268-2, BCAT2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68586-1-PBS. Detection antibody: 68586-4-PBS. Standard: Ag9595. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP50268-3, BCAT2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68586-5-PBS. Detection antibody: 68586-4-PBS. Standard: Ag9595. Range: 0.781-100 ng/mL.