

For Research Use Only

BCKDHB Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 68596-4-PBS



Basic Information

Catalog Number:

68596-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG3

Immunogen Catalog Number:

AG33802

GenBank Accession Number:

BC040139

GeneID (NCBI):

594

UNIPROT ID:

P21953

Full Name:

branched chain keto acid
dehydrogenase E1, beta polypeptide

Calculated MW:

392 aa, 43 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

2C5C12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50570-2, BCKDHB Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68596-3-PBS. Detection antibody: 68596-4-PBS. Standard: Ag33802. Range: 0.781-50 ng/mL



Cytometric bead array standard curve of MP50570-3, BCKDHB Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68596-4-PBS. Detection antibody: 68596-1-PBS. Standard: Ag33802. Range: 1.563-100 ng/mL