

For Research Use Only

# SLC1A4 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68601-3-PBS

## Basic Information

Catalog Number:

68601-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG17328

GenBank Accession Number:

BC026216

GeneID (NCBI):

6509

UNIPROT ID:

P43007

Full Name:

solute carrier family 1  
(glutamate/neutral amino acid  
transporter), member 4

Calculated MW:

532 aa, 56 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

2E3B8

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50199-1, SLC1A4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68601-2-PBS. Detection antibody: 68601-3-PBS. Standard: Ag17328. Range: 0.098-12.5 ng/mL