

For Research Use Only

# ACSS1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68608-4-PBS



## Basic Information

<b>Catalog Number:</b> 68608-4-PBS	<b>GenBank Accession Number:</b> BC039261	<b>Purification Method:</b> Protein G Magarose purification
<b>Size:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 84532	<b>CloneNo.:</b> 4A6C10
<b>Source:</b> Mouse	<b>UNIPROT ID:</b> Q9NUB1	
<b>Isotype:</b> IgG1	<b>Full Name:</b> acyl-CoA synthetase short-chain family member 1	
<b>Immunogen Catalog Number:</b> AG10853	<b>Calculated MW:</b> 689 aa, 75 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50764-2, ACS1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68608-2-PBS. Detection antibody: 68608-4-PBS. Standard: Ag10853. Range: 1.563-100 ng/mL.