

For Research Use Only

RNF160 Monoclonal antibody, PBS Only



Catalog Number: 68610-1-PBS

Basic Information

Catalog Number: 68610-1-PBS	GenBank Accession Number: BC140790	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 26046	CloneNo.: 1D2A2
Source: Mouse	UNIPROT ID: O94822	
Isotype: IgG1	Full Name: ring finger protein 160	
Immunogen Catalog Number: AG29471	Observed MW: 180-200 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human, rat, mouse

Background Information

The E3 ubiquitin ligase, Ltn1 (also named RNF160), is a key regulator of RQC that targets substrates for ubiquitination (PMID: 35981909).

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68610-1-Ig (RNF160 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68610-1-PBS in a different storage buffer formulation.