

For Research Use Only

PRODH Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68621-5-PBS



Basic Information

Catalog Number:

68621-5-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG18429

GenBank Accession Number:

BC118597

GeneID (NCBI):

5625

UNIPROT ID:

O43272

Full Name:

proline dehydrogenase (oxidase) 1

Calculated MW:

600 aa, 68 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1G3C5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50609-2, PRODH Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68621-4-PBS. Detection antibody: 68621-5-PBS. Standard:Ag18429. Range: 0.391-100 ng/mL.