

For Research Use Only

PRODH Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68621-6-PBS



Basic Information

Catalog Number: 68621-6-PBS	GenBank Accession Number: BC118597	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 5625	CloneNo.: 2E11G6
Source: Mouse	UNIPROT ID: O43272	
Isotype: IgG2a	Full Name: proline dehydrogenase (oxidase) 1	
Immunogen Catalog Number: AG18429	Calculated MW: 600 aa, 68 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50609-3, PRODH Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68621-6-PBS. Detection antibody: 68621-7-PBS. Standard: Ag18429. Range: 0.391-6.25 ng/mL