

For Research Use Only

Cytokeratin 8 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68627-3-PBS



Basic Information

Catalog Number:

68627-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG0488

GenBank Accession Number:

BC000654

GeneID (NCBI):

3856

UNIPROT ID:

P05787

Full Name:

keratin 8

Calculated MW:

54 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1E10E4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51080-2, Cytokeratin 8 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68627-3-PBS. Detection antibody: 68627-4-PBS. Standard: Ag0488. Range: 0.391-100 ng/mL