

For Research Use Only

# TFEB Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68632-2-PBS

## Basic Information

Catalog Number: 68632-2-PBS	GenBank Accession Number: BC032448	Purification Method: Protein G purification
Size: 1mg/ml	GeneID (NCBI): 7942	CloneNo.: 1D1B9
Source: Mouse	UNIPROT ID: P19484	
Isotype: IgG1	Full Name: transcription factor EB	
Immunogen Catalog Number: AG24850	Calculated MW: 476 aa, 53 kDa	

## Applications

Tested Applications:  
Indirect ELISA, Cytometric bead array

Species Specificity:  
Human

## Background Information

### Storage

Storage:  
Store at -80°C.  
Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50099-1, TFEB Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68632-1-PBS. Detection antibody: 68632-2-PBS. Standard:Ag24850. Range: 0.195-50 ng/mL