

For Research Use Only

# PGF Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68644-1-PBS

## Basic Information

Catalog Number:

68644-1-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG2b

GenBank Accession Number:

NP\_001193941.1

GeneID (NCBI):

5228

UNIPROT ID:

P49763-2

Full Name:

placental growth factor

Calculated MW:

17KD

Purification Method:

Protein A purification

CloneNo.:

3E9D2

Affinity:

$K_D = 2.49 \times 10^{-9}$

$K_{off} = 1.28 \times 10^{-4}$

$K_{on} = 5.12 \times 10^4$

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sandwich ELISA standard curve of MP50063-1, PGF Monoclonal Matched Antibody Pair - PBS only. 68644-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg0653. 68644-1-PBS was HRP conjugated as the detection antibody. Range: 31.25-2000 pg/mL



Sandwich ELISA standard curve of MP50063-2, PGF Monoclonal Matched Antibody Pair - PBS only. 68644-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg0653. 68644-1-PBS was HRP conjugated as the detection antibody. Range: 31.25-2000 pg/mL