

For Research Use Only

# DPF3 Monoclonal antibody, PBS Only



Catalog Number: 68646-1-PBS

## Basic Information

Catalog Number:

68646-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG33877

GenBank Accession Number:

BC060801

GeneID (NCBI):

8110

UNIPROT ID:

Q92784

Full Name:

D4, zinc and double PHD fingers,  
family 3

Calculated MW:

357 aa, 40 kDa

Observed MW:

43 kDa

Purification Method:

Protein A purification

CloneNo.:

1F3H4

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human, rabbit, rat, mouse, pig

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68646-1-Ig (DPF3 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68646-1-PBS in a different storage buffer formulation.