

For Research Use Only

GOLGA1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68648-3-PBS



Basic Information

Catalog Number:

68648-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG33790

GenBank Accession Number:

BC032853

GeneID (NCBI):

2800

UNIPROT ID:

Q92805

Full Name:

golgi autoantigen, golgin subfamily a, 1

Calculated MW:

767 aa, 88 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1H2C5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50871-2, GOLGA1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68648-3-PBS. Detection antibody: 68648-4-PBS. Standard:Ag33790. Range: 1.563-100 ng/mL.



Cytometric bead array standard curve of MP50871-3, GOLGA1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68648-3-PBS. Detection antibody: 68648-1-PBS. Standard:Ag33790. Range: 1.563-100 ng/mL.