

For Research Use Only

RS1 Monoclonal antibody

Catalog Number: 68654-1-Ig



Basic Information

| | | |
|---|--|---|
| Catalog Number: 68654-1-Ig | GenBank Accession Number: BC141638 | Purification Method: Protein G purification |
| Size: 1000 µg/ml | GeneID (NCBI): 6247 | CloneNo.: 1A2D4 |
| Source: Mouse | UNIPROT ID: O15537 | Recommended Dilutions: WB 1:5000-1:50000 |
| Isotype: IgG1 | Full Name: retinoschisin 1 | |
| Immunogen Catalog Number: AG18885 | Calculated MW: 26 kDa | |
| | Observed MW: 26 kDa | |

Applications

| | |
|--|--|
| Tested Applications: WB, ELISA | Positive Controls: WB : rabbit retina tissue, mouse retina tissue, rat retina tissue |
| Species Specificity: Human, rat, mouse, rabbit | |

Background Information

RS1 is also named as XLRS1. RS1 can bind negatively charged membrane lipids, such as phosphatidylserine and phosphoinositides. RS1 is required for normal structure and function of the retina (PMID:19093009). It may play a role in cell-cell adhesion processes in the retina (PMID:27114531). RS1 is high expression in the retina (at protein level) (PMID:10915776, PMID:9326935). It is detected in the inner segment of the photoreceptors, the inner nuclear layer, the inner plexiform layer and the ganglion cell layer. At the macula, it is expressed in both the outer and inner nuclear layers and in the inner plexiform layer (PMID:10915776, PMID:10915776). And RS1 is undetectable in the inner plexiform layers and the inner nuclear layer (PMID:10915776)

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68654-1-Ig (RS1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.