

For Research Use Only

# IL-6 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68659-3-PBS

## Basic Information

Catalog Number:

68659-3-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

HZ-1019

GenBank Accession Number:

GeneID (NCBI):

3569

ENSEMBL Gene ID:

ENSG00000136244

Full Name:

interleukin 6 (interferon, beta 2)

Purification Method:

Protein G purification

CloneNo.:

4D12G2

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50032-2, IL-6 Monoclonal Matched Antibody Pair - BSA and Azide free. Capture antibody: 68659-3-PBS. Detection antibody: 66146-2-PBS. Standard:HZ-1019 . Range: 1.563-100 ng/mL.