

For Research Use Only

# HGF Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68660-2-PBS

## Basic Information

Catalog Number:

68660-2-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

HZ-1084

GenBank Accession Number:

GeneID (NCBI):

3082

Full Name:

hepatocyte growth factor  
(hepapoietin A; scatter factor)

Purification Method:

Protein G purification

CloneNo.:

2E5G2

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sandwich ELISA standard curve of MP50034-1, HGF Monoclonal Matched Antibody Pair - PBS only. 68660-1-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard HZ-1084. 68660-2-PBS was HRP conjugated as the detection antibody. Range: 0.313-20 ng/mL