

For Research Use Only

FLT3 Ligand Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68664-1-PBS

Basic Information

Catalog Number:

68664-1-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

HZ-1151

GenBank Accession Number:

GeneID (NCBI):

2323

Full Name:

fms-related tyrosine kinase 3 ligand

Purification Method:

Protein A purification

CloneNo.:

1C1H5

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50046-1, FLT3LG Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68664-1-PBS. Detection antibody: 68664-2-PBS. Standard: HZ-1151. Range: 15.6-1000 pg/mL.