

For Research Use Only

IL-28B Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68668-1-PBS

Basic Information

Catalog Number:

68668-1-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

HZ-1245

GenBank Accession Number:

GeneID (NCBI):

282617

Full Name:

interleukin 28B (interferon, lambda 3)

Purification Method:

Protein G purification

CloneNo.:

3G6E3

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50007-1, IL-28B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68668-1-PBS. Detection antibody: 68668-2-PBS. Standard: HZ-1245. Range: 15.6-1000 pg/mL.