

For Research Use Only

PCSK9 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68673-2-PBS

Basic Information

Catalog Number: 68673-2-PBS	GenBank Accession Number: NM_174936.4	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 255738	CloneNo.: 3H8H9
Source: Mouse	UNIPROT ID: Q8NBP7	
Isotype: IgG2a	Full Name: proprotein convertase subtilisin/kexin type 9	
	Calculated MW: 74 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50055-1, PCSK9 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68673-1-PBS. Detection antibody: 68673-2-PBS. Standard: Eg0199. Range: 31.25-2000 pg/mL