

For Research Use Only

VISTA Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68676-1-PBS



Basic Information

Catalog Number: 68676-1-PBS	GenBank Accession Number: BC020568	Purification Method: Protein G purification
Size: 1mg/ml	GeneID (NCBI): 64115	CloneNo.: 2B3H6
Source: Mouse	UNIPROT ID: Q9H7M9	
Isotype: IgG1	Full Name: chromosome 10 open reading frame 54	
	Calculated MW: 311 aa, 34 kDa	

Applications

Tested Applications:
Cytometric bead array, Sandwich ELISA, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50058-1, VISTA Monoclonal Matched Antibody Pair - PBS only. 68676-1-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg0017. 68676-2-PBS was HRP conjugated as the detection antibody. Range: 0.156-10 ng/mL