

For Research Use Only

FCGR3B Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68677-1-PBS

Basic Information

Catalog Number:

68677-1-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG2a

GenBank Accession Number:

NM_000570

GeneID (NCBI):

2215

UNIPROT ID:

O75015

Full Name:

Fc fragment of IgG, low affinity IIIb, receptor (CD16b)

Calculated MW:

26 kDa

Purification Method:

Protein A purification

CloneNo.:

3C4G3

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50059-1, FCGR3B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68677-1-PBS. Detection antibody: 68677-2-PBS. Standard: Eg0308. Range: 0.39-25 ng/mL.