

For Research Use Only

FCGR3B Monoclonal antibody, PBS Only (Capture/Detector)



Catalog Number: 68677-2-PBS

Basic Information

Catalog Number: 68677-2-PBS	GenBank Accession Number: NM_000570	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 2215	CloneNo.: 4E7D3
Source: Mouse	UNIPROT ID: O75015	
Isotype: IgG2a	Full Name: Fc fragment of IgG, low affinity IIIb, receptor (CD16b)	
	Calculated MW: 26 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50059-1, FCGR3B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68677-1-PBS. Detection antibody: 68677-2-PBS. Standard: Eg0308. Range: 0.39-25 ng/mL



Cytometric bead array standard curve of MP50059-2, FCGR3B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68677-2-PBS. Detection antibody: 68677-3-PBS. Standard: Eg0308. Range: 0.195-100 ng/mL