

For Research Use Only

# IL1RL1 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68680-2-PBS

## Basic Information

Catalog Number:

68680-2-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

BC030975

GeneID (NCBI):

9173

UNIPROT ID:

Q01638-2

Full Name:

interleukin 1 receptor-like 1

Calculated MW:

63 kDa, 37 kDa, 30 kDa

Purification Method:

Protein G purification

CloneNo.:

2F2C3

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sandwich ELISA standard curve of MP50068-1, IL1RL1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68680-1-PBS. Detection antibody: 68680-2-PBS. Standard: Eg0116. Range: 31.25-2000 pg/mL