

For Research Use Only

TNFR1 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68685-2-PBS

Basic Information

Catalog Number:

68685-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

BC010140

GeneID (NCBI):

7132

UNIPROT ID:

P19438

Full Name:

tumor necrosis factor receptor
superfamily, member 1A

Calculated MW:

455 aa, 50 kDa

Purification Method:

Protein G purification

CloneNo.:

1D12C3

Applications

Tested Applications:

Cytometric bead array, ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Standard curve of MP50042-1, TNFR1 Monoclonal Matched Antibody Pair - BSA and Azide free, measured by sandwich ELISA. Capture antibody: 68685-1-PBS. Detection antibody: 68685-2-PBS (HRP conjugated). Standard: Eg0531. Range: 7.8-500 pg/mL