

For Research Use Only

CLU Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68687-1-PBS

Basic Information

Catalog Number:

68687-1-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

BC010514

GeneID (NCBI):

1191

UNIPROT ID:

P10909

Full Name:

clusterin

Calculated MW:

449 aa, 53 kDa

Purification Method:

Protein G purification

CloneNo.:

2G4A9

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50075-1, CLU
Monoclonal Matched Antibody Pair, PBS Only.
Capture antibody: 68687-1-PBS. Detection
antibody: 68687-2-PBS. Standard: Eg0225. Range:
0.313-20 ng/mL