

For Research Use Only

S100A14 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68690-1-PBS

Basic Information

Catalog Number:

68690-1-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG29970

GenBank Accession Number:

BC005019

GeneID (NCBI):

57402

UNIPROT ID:

Q9HCY8

Full Name:

S100 calcium binding protein A14

Calculated MW:

12 kDa

Purification Method:

Protein G purification

CloneNo.:

1D12A2

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50082-1, S100A14 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68690-1-PBS. Detection antibody: 68690-2-PBS. Standard: Ag29970. Range: 0.156-10 ng/mL