

For Research Use Only

# TIMP1 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68691-4-PBS

## Basic Information

<b>Catalog Number:</b> 68691-4-PBS	<b>GenBank Accession Number:</b> BC000866	<b>Purification Method:</b> Protein G purification
<b>Size:</b> 1mg/ml	<b>GeneID (NCBI):</b> 7076	<b>CloneNo.:</b> 2F2C2
<b>Source:</b> Mouse	<b>UNIPROT ID:</b> P01033	<b>Affinity:</b> $K_D = 1.57 \times 10^{-9}$
<b>Isotype:</b> IgG1	<b>Full Name:</b> TIMP metalloproteinase inhibitor 1	$K_{off} = 3.48 \times 10^{-4}$ $K_{on} = 2.22 \times 10^5$
	<b>Calculated MW:</b> 23 kDa	

## Applications

**Tested Applications:**  
Sandwich ELISA, Indirect ELISA

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sandwich ELISA standard curve of MP50084-2, TIMP1 Monoclonal Matched Antibody Pair - PBS only. 68691-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg0638. 68691-4-PBS was HRP conjugated as the detection antibody. Range: 62.5-4000 pg/mL.