

For Research Use Only

# gp130/IL6ST Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68692-2-PBS

## Basic Information

Catalog Number:

68692-2-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

NM\_002184.3

GeneID (NCBI):

3572

UNIPROT ID:

P40189-1

Full Name:

interleukin 6 signal transducer  
(gp130, oncostatin M receptor)

Calculated MW:

104 kDa

Purification Method:

Protein G purification

CloneNo.:

1G9E1

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Sandwich ELISA standard curve of MP50086-1, gp130/IL6ST Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68692-1-PBS. Detection antibody: 68692-2-PBS. Standard: Eg0654. Range: 31.25-2000 pg/mL