

For Research Use Only

ICOSLG Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68694-2-PBS

Basic Information

Catalog Number:

68694-2-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

BC064637

GeneID (NCBI):

23308

UNIPROT ID:

O75144

Full Name:

inducible T-cell co-stimulator ligand

Calculated MW:

33 kDa, 34 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3A8C12

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50083-1, ICOSLG Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68694-1-PBS. Detection antibody: 68694-2-PBS. Standard: Eg0008. Range: 15.6-1000 pg/mL.