

For Research Use Only

Rat Mmp9 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68695-3-PBS

Basic Information

Catalog Number:

68695-3-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

NM_031055.2

GeneID (NCBI):

81687

Full Name:

matrix metalloproteinase 9

Calculated MW:

79kd

Purification Method:

Protein G purification

CloneNo.:

2C6B6

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

Rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50079-1, Rat Mmp9 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68695-1-PBS. Detection antibody: 68695-3-PBS. Standard: Eg0439. Range: 62.5-4000 pg/mL.