

For Research Use Only

TF Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68697-2-PBS

Basic Information

Catalog Number:

68697-2-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

BC059367

GeneID (NCBI):

7018

UNIPROT ID:

P02787

Full Name:

transferrin

Calculated MW:

698 aa, 77 kDa

Purification Method:

Protein G purification

CloneNo.:

1G5E1

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50076-1, TF Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68697-1-PBS. Detection antibody: 68697-2-PBS. Standard: Eg0114. Range: 23.4-1500 pg/mL.