

For Research Use Only

IL2RA Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68698-1-PBS

Basic Information

Catalog Number:

68698-1-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

BC156627

GeneID (NCBI):

3559

UNIPROT ID:

P01589

Full Name:

interleukin 2 receptor, alpha

Calculated MW:

272 aa, 31 kDa

Purification Method:

Protein G purification

CloneNo.:

1A1B2

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50073-1, IL2RA Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68698-1-PBS. Detection antibody: 68698-2-PBS. Standard: Eg0408. Range: 78.1-2500 pg/mL.