

For Research Use Only

PDGF-D Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68701-1-PBS



Basic Information

| | | |
|---|---|--|
| Catalog Number: 68701-1-PBS | GenBank Accession Number: BC030645 | Purification Method: Protein G Magarose purification |
| Size: 1 mg/ml | GeneID (NCBI): 80310 | CloneNo.: 2A4B7 |
| Source: Mouse | UNIPROT ID: Q9GZP0 | |
| Isotype: IgG1 | Full Name: platelet derived growth factor D | |
| Immunogen Catalog Number: AG29199 | Calculated MW: 43 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50964-1, PDGFD Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68701-2-PBS. Detection antibody: 68701-1-PBS. Standard:Ag5154. Range: 0.098-100 ng/mL



Cytometric bead array standard curve of MP50964-2, PDGFD Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68701-3-PBS. Detection antibody: 68701-1-PBS. Standard:Ag5154. Range: 0.098-100 ng/mL