

For Research Use Only

LEPR Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68703-1-PBS

Basic Information

Catalog Number: 68703-1-PBS	GenBank Accession Number: NM_002303	Purification Method: Protein G purification
Size: 1mg/ml	GeneID (NCBI): 3953	CloneNo.: 3A4E5
Source: Mouse	UNIPROT ID: P48357	
Isotype: IgG1	Full Name: leptin receptor	
Immunogen Catalog Number: AG29621	Calculated MW: 132 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50066-1, LEPR
Monoclonal Matched Antibody Pair, PBS Only.
Capture antibody: 68703-1-PBS. Detection
antibody: 68703-2-PBS. Standard: Ag29621. Range:
39.06-2500 pg/mL