

For Research Use Only

# AGT Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68704-1-PBS

## Basic Information

<b>Catalog Number:</b> 68704-1-PBS	<b>GenBank Accession Number:</b> BC011519	<b>Purification Method:</b> Protein G purification
<b>Size:</b> 1mg/ml	<b>GeneID (NCBI):</b> 183	<b>CloneNo.:</b> 2G2G7
<b>Source:</b> Mouse	<b>UNIPROT ID:</b> P01019	
<b>Isotype:</b> IgG1	<b>Full Name:</b> angiotensinogen (serpin peptidase inhibitor, clade A, member 8)	
<b>Immunogen Catalog Number:</b> AG16658	<b>Calculated MW:</b> 485 aa, 53 kDa	

## Applications

**Tested Applications:**  
Sandwich ELISA, Indirect ELISA

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sandwich ELISA standard curve of MP50065-1, AGT Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68704-1-PBS. Detection antibody: 68704-2-PBS. Standard: Ag16658. Range: 0.156-10 ng/mL