

For Research Use Only

NPTX2 Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 68709-3-PBS



Basic Information

Catalog Number: 68709-3-PBS	GenBank Accession Number: BC009924	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 4885	CloneNo.: 3A12G3
Source: Mouse	UNIPROT ID: P47972	
Isotype: IgG1	Full Name: neuronal pentraxin II	
Immunogen Catalog Number: AG33808	Calculated MW: 47 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50983-1, NPTX2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68709-2-PBS. Detection antibody: 68709-3-PBS. Standard:Ag33808. Range: 0.098-100 ng/mL.



Cytometric bead array standard curve of MP50983-2, NPTX2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68709-3-PBS. Detection antibody: 68709-4-PBS. Standard:Ag33808. Range: 0.391-100 ng/mL.